


■ companion ■

Compatibility – Companion power add-on for manual wheelchairs



COMPATIBLE WITH YOUR


Rigid Manual Wheelchair

The Companion is compatible with most rigid manual wheelchairs that have a straight camber tube and crossbar diameter between 19mm (0.75") and 50.8mm (2").

COMPATIBLE WITH YOUR

Folding Manual Wheelchair

The Companion is compatible with most foldable manual wheelchairs that have a frame tube diameter between 19mm (0.75") and 50.8mm (2").



If your wheelchair has a different bar or frame tube diameter, please [contact us](#) for a customized solution.

Exceed your reach with your Companion power add-on.