

companion



Technical Specs - Companion power add-on for manual wheelchairs

Motor	250W Hub motor with 8 inch solid tire
Brakes	Powerful drum brakes
Battery	Lithium battery 36V 7.5A
Range	Up to 25km/16miles on full charge*
Power Switch	ON/OFF Power switch with two keys
Turning Radius	32 inches
Handlebar angle range	180 degrees
Speed	Three speed mode control - maximum 12km/h
Display	2inch LCD display - battery charge status, trip distance and speed
Weight - Companion Unit	9.4kgs/20lbs**
Weight - Adaptor/Coupler	0.8kgs/1.8lbs

^{*}Range measured in laboratory conditions: flat clean surface, no inclines, 150lbs dummy. Actual range might vary depending on road, incline, and weather conditions.

Exceed your reach with your Companion power add-on.

^{**}For Companion optimal performance, total combined wheelchair and user weight should not exceed 280lbs/130kgs.