




Front wrap options



## Technical Specs – Companion power assist

- Motor -250W Hub motor with 8 inch solid tire
- Brakes – Powerful drum brakes
- Battery– Lithium battery 36V 7.5A
- Range - Up to 25km/16miles on full charge
- Power Switch – ON/OFF Power switch with 2 keys
- Turning Radius – 32 inches
- Handlebar angle range – 180 degrees
- Speed - Three speed mode control - max 12km/h
- Display–2inch LCD display - battery charge status, trip distance and speed
- Weight - Companion Unit–9.4kgs/20lbs
- Weight - Adaptor/Coupler - 0.8kgs/1.8lbs

Exceed your reach with your Companion power assist.



### CHEELCARE

16 Sims Crescent, unit 20  
Richmond Hill, Ontario  
Canada, L4B 2P1

Phone: 647-800 2680  
E-mail: info@cheelcare.com  
www.cheelcare.com

## companion

Power Assist for Manual Wheelchairs



EXCEED YOUR REACH

**COMPATIBLE** WITH YOUR  
Rigid Manual  
Wheelchair



**COMPATIBLE** WITH YOUR  
Folding Manual  
Wheelchair



**COMPATIBLE** WITH YOUR  
Pediatric Manual  
Wheelchair



- Lightweight – less than 10 kgs – Ideal portable solution for car, sea and air travel
- Outstanding maneuverability – Turns in place and easily steers through narrow spaces
- Sleek Compact design based on 8" wheel
- Unique Quick-Release connect and disconnect mechanism – less than 10 seconds!
- Smooth, powerful and fast – 3 speeds – comfortable ride on any occasion and place
- Long range – 25 kms and fast 3-hour charge
- Safe – Amazingly stable, safe and secure. Powerful drum brakes.
- Connect to most RIGID and FOLDING manual wheelchairs, both adult and pediatric. Endless adjustment options.
- Highly customizable, make it look Yours! (*shown in front wrapping options*)
- Your perfect Companion all day, every day. At school, at work, at the mall; visiting friends or your doctor. Always with you.



Exceed your reach with your  
Companion power assist.

**Companion Accessories**

Detachable  
basket



Throttle Assist

Q-type bell



Lights



Brakes with  
park assist



**cheelcare**<sup>TM</sup>  
exceed your reach

EXCEED YOUR REACH